

Product Information

Control of Non-PEGylated Liposome with Doxorubicin (2.0mg, mL)

Cat. No.: **MPL0944K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of Hydrogenated Soy PC and Cholesterol (60:40). The mean particle size is 100 nm and the lipid concentration is 24.61mM. It should be used as a negative control for Non-PEGylated Liposomal with Doxorubicin (2.0mg/mL).

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL; 5 mL

Storage

2-8°C

Storage Shelf Time

4 months