

# Product Information

## Succinyl Liposome with Doxorubicin (2.0mg, mL)

Cat. No.: **MPL0946K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of Hydrogenated Soy PC, Cholesterol, DSPE-PEG(2000) and DSPE-PEG(2000)-Succinyl (57:38:4:1). The mean particle size is 100 nm and the lipid concentration is 21.58mM.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL; 5 mL

### Storage

2-8°C

### Storage Shelf Time

4 months