

Product Information

Cyanur Liposome with Doxorubicin (2.0mg, mL)

Cat. No.: **MPL0947K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of Hydrogenated Soy PC, Cholesterol, DSPE-PEG(2000) and DSPE-PEG(2000)-Cyanur (57:38:4:1). The mean particle size is 100 nm and the lipid concentration is 21.58mM.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL; 5 mL

Storage

2-8°C

Storage Shelf Time

4 months