

# Product Information

## NTA(Ni) Liposome

Cat. No.: **MPL0971K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of L-alpha-Phosphatidylcholine, Cholesterol and 1,2-Dioleoyl-sn-glycero-3-[(N-(5-amino-1-carboxypentyl)iminodiacetic acid)succinyl] (nickel salt) (69:30:1). The mean particle size is 100 nm and the lipid concentration is 22.45mM.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL; 5 mL

### Storage

2-8°C

### Storage Shelf Time

4 months