

Product Information

Carboxylic Acid PEGylated Fluorescent Liposome

Cat. No.: **MPL0979K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of Hydrogenated Soy PC, Cholesterol, DSPE-PEG(2000) and DSPE-PEG(2000)-Carboxylic Acid (57:38:4:1). The mean particle size is 100 nm and the lipid concentration is 21.55mM. Fluorescent dyes can be selected based on the requirements.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

2-8°C

Storage Shelf Time

4 months