

# Product Information

## Biotin Fluorescent Liposome

Cat. No.: **MPL0981K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of L-alpha-Phosphatidylcholine, Cholesterol and 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-(biotinyl) (sodium salt) (69:30:1). The mean particle size is 100 nm and the lipid concentration is 22.45mM. Fluorescent dyes can be selected based on the requirements.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL

### Storage

2-8°C

### Storage Shelf Time

4 months