

Product Information

Azide Fluorescent Liposome

Cat. No.: **MPL0986K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of L-alpha-Phosphatidylcholine, Cholesterol and 1,2-Dipalmitoyl-sn-glycero-3-phosphoethanolamine-N-(6-azidohehexanoyl) (69:30:1). The mean particle size is 100 nm and the lipid concentration is 22.45mM. Fluorescent dyes can be selected based on the requirements.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

2-8°C

Storage Shelf Time

4 months