

Product Information

PG-based DBCO Fluorescent Liposome

Cat. No.: **MPL1001K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 1,2-dioleoyl-sn-glycero-3-phospho-(1'-rac-glycerol) (sodium salt), 1,2-dioleoyl-sn-glycero-3-phosphocholine, Cholesterol and 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-dibenzocyclooctyl (20:49:30:1). The mean particle size is 100 nm and the lipid concentration is 10mM. Fluorescent dyes can be selected based on the requirements.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months