

Product Information

DOPC, DOPG Liposomes Co-encapsulating NBD, Rhodamine

Cat. No.: **MPL1009K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of DOPC and DOPG (DOPC:DOPG=90:10, 80:20, 70:30, 60:40, 50:50, 40:60, 30:70, 20:80, 10:90, 0:100). The mean particle size is 100 nm. Fluorescent dyes can be selected based on the requirements.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months