

Product Information

POPC, Chol, POPG Liposomes Co-encapsulating NBD, Rhodamine

Cat. No.: **MPL1010K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of POPC, Chol and POPG (POPC:Chol:POPG=10:4:1). The mean particle size is 100 nm. Fluorescent dyes can be selected based on the requirements.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months