

# Product Information

## POPE, POPS, POPC Liposomes Co-encapsulating NBD, Rhodamine

Cat. No.: **MPL1013K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of POPE, POPS and POPC (POPE:POPS:POPC=5:5:2). The mean particle size is 100 nm. Fluorescent dyes can be selected based on the requirements.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL

### Storage

4°C, avoiding light

### Storage Shelf Time

2 months