

Product Information

DOTAP Liposomes Co-encapsulating NBD, Rhodamine

Cat. No.: **MPL1015K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 100% DOTAP or DOTAP:Chol (1:1) or DOTAP:DOPE (1:1). The mean particle size is 100 nm. Fluorescent dyes can be selected based on the requirements.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months