

# Product Information

## POPC, Chol, POPG Liposomes Containing Octadecyl Rhodamine B

Cat. No.: **MPL1016K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 1-palmitoyl-2-oleoyl-sn-glycero-3-phosphocholine, Cholesterol and 1-palmitoyl-2-oleoyl-sn-glycero-3-phospho-(1'-rac-glycerol) (sodium salt) (10:4:1). The mean particle size is 100 nm and the lipid concentration is 10mM. Fluorescent dye is Octadecyl Rhodamine B Chloride.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL

### Storage

4°C, avoiding light

### Storage Shelf Time

2 months