

Product Information

POPE, POPG Liposomes Containing Octadecyl Rhodamine B

Cat. No.: **MPL1018K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 1-palmitoyl-2-oleoyl-sn-glycero-3-phosphoethanolamine and 1-palmitoyl-2-oleoyl-sn-glycero-3-phospho-(1'-rac-glycerol) (sodium salt) (3:1). The mean particle size is 100 nm and the lipid concentration is 10mM. Fluorescent dye is Octadecyl Rhodamine B Chloride.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months