

Product Information

DOPC, DOPS Liposomes Containing Octadecyl Rhodamine B

Cat. No.: **MPL1019K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of DOPC and DOPS (DOPC:DOPS=90:10, 80:20, 70:30, 60:40, 50:50, 40:60, 30:70, 20:80, 10:90, 0:100). The mean particle size is 100 nm. Fluorescent dye is Octadecyl Rhodamine B Chloride.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months