

Product Information

DOTAP Liposomes Containing Octadecyl Rhodamine B

Cat. No.: **MPL1022K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 100% DOTAP or DOTAP:Chol (1:1) or DOTAP:DOPE (1:1). The mean particle size is 100 nm. Fluorescent dye is Octadecyl Rhodamine B Chloride.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months