

# Product Information

## DOPC, DOPG Liposomes Containing Calcein or Carboxyfluorescein

Cat. No.: **MPL1025K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of DOPC and DOPG (DOPC:DOPG=90:10, 80:20, 70:30, 60:40, 50:50, 40:60, 30:70, 20:80, 10:90, 0:100). The mean particle size is 100 nm. Fluorescent dye is 5(6)-Carboxyfluorescein or Calcein.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL

### Storage

4°C, avoiding light

### Storage Shelf Time

2 months