

Product Information

POPE, POPS, POPC Liposomes Containing Calcein or Carboxyfluorescein

Cat. No.: **MPL1026K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 1-palmitoyl-2-oleoyl-sn-glycero-3-phosphoethanolamine, 1-palmitoyl-2-oleoyl-sn-glycero-3-phospho-L-serine (sodium salt) and 1-palmitoyl-2-oleoyl-sn-glycero-3-phosphocholine (5:5:2). Fluorescent dye is 5(6)-Carboxyfluorescein or Calcein.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL

Storage

4°C, avoiding light

Storage Shelf Time

2 months