

# Product Information

## DOPG, DOPC Plain Liposomes (0.5:99.5)

Cat. No.: **MPL1215K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of DOPG, DOPC. The available particle size is 50 nm, 100 nm, 200 nm, 1-2  $\mu$ m.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

5 ml

### Storage

4°C, avoiding light

### Storage Shelf Time

4 months