

Product Information

18:2 (Cis) PC (DLPC) Plain Liposomes

Cat. No.: **MPL1307K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 100% 18:2 (Cis) PC (DLPC). The available particle size is 50 nm, 100 nm, 200 nm, 1-2 µm.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

5 ml

Storage

4°C, avoiding light

Storage Shelf Time

4 months