

# Product Information

## 18:0-20:4 PC, Cholesterol Plain Liposomes

Cat. No.: **MPL1336K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of 18:0-20:4 PC, Cholesterol (70:30 ). The available particle size is 50 nm, 100 nm, 200 nm, 1-2  $\mu$ m.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

5 ml

### Storage

4°C, avoiding light

### Storage Shelf Time

4 months