

Product Information

Human MSP1D1 (D73C), His

Cat. No.: **MPN0033K**

This product is for research use only and is not intended for diagnostic use.

This product is a MSP protein for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Human

Protein Tag

His tag

Form

Lyophilized powder

Biochem Actions

9-10 nm in diameter Nanodiscs

Size

2 mg; 10 mg

Storage

Aliquot and store at -20°C or lower.