

# Product Information

## Human MSP1D1 $\Delta$ H5-His with DMPC

Cat. No.: **MPN0038K**

This product is for research use only and is not intended for diagnostic use.

This product is a 0.5 mM concentrated pre-assembled nanodisc for cell-free protein expression. It is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Type

Nanodiscs

### Species

Human

### Phospholipids

DMPC

### Protein Tag

His tag

### Form

Assembled

### Biochem Actions

7-8 nm in diameter Nanodiscs

### Size

50  $\mu$ L

### Storage

Aliquot and store at -20°C or lower.