

Product Information

Human MSP1D1 Δ H5-His with POPC

Cat. No.: **MPN0039K**

This product is for research use only and is not intended for diagnostic use.

This product is a 0.5 mM concentrated pre-assembled nanodisc for cell-free protein expression. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Human

Phospholipids

POPC

Protein Tag

His tag

Form

Assembled

Biochem Actions

7-8 nm in diameter Nanodiscs

Size

50 μ L

Storage

Aliquot and store at -20°C or lower.