

Product Information

Human MSP1D1 Δ H5-tag-free with DMPG

Cat. No.: **MPN0046K**

This product is for research use only and is not intended for diagnostic use.

This product is a 0.5 mM concentrated pre-assembled nanodisc for cell-free protein expression. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Human

Phospholipids

DMPG

Protein Tag

Tag-free

Form

Assembled

Biochem Actions

7-8 nm in diameter Nanodiscs

Size

50 μ L

Storage

Aliquot and store at -20°C or lower.