

Product Information

Human MSP1D1 Δ H5-His with DMPC Biotinyl PE

Cat. No.: **MPN0061K**

This product is for research use only and is not intended for diagnostic use.

This product is a 0.5 mM concentrated pre-assembled nanodisc for cell-free protein expression. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Human

Phospholipids

Biotin-PE, DMPC

Protein Tag

His tag

Form

Assembled

Biochem Actions

7-8 nm in diameter Nanodiscs

Size

50 μ L

Storage

Aliquot and store at -20°C or lower.