

Product Information

Human MSP1E3D1-His with DMPC Biotinyl PE

Cat. No.: **MPN0063K**

This product is for research use only and is not intended for diagnostic use.

This product is a 0.5 mM concentrated pre-assembled nanodisc for cell-free protein expression. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Human

Phospholipids

Biotin-PE, DMPC

Protein Tag

His tag

Form

Assembled

Biochem Actions

12-14 nm in diameter Nanodiscs

Size

50 µL

Storage

Aliquot and store at -20°C or lower.