

Product Information

Mouse MSP1E3D1-His with DMPC

Cat. No.: **MPN0068K**

This product is for research use only and is not intended for diagnostic use.

This product is a 0.5 mM concentrated pre-assembled nanodisc for cell-free protein expression. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Mouse

Phospholipids

DMPC

Protein Tag

His tag

Form

Assembled

Biochem Actions

12-14 nm in diameter Nanodiscs

Size

50 µL

Storage

Aliquot and store at -20°C or lower.