

Product Information

Human MSP1E3D1-His with DMPC Nanodisc Assembly Kit

Cat. No.: **MPN0089K**

This product is for research use only and is not intended for diagnostic use.

This product is a nanodisc assembly kit to prepare empty nanodiscs, including MSP protein, phospholipids and sodium cholate. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Nanodiscs

Species

Human

Phospholipids

DMPC

Protein Tag

His tag

Form

Kit

Biochem Actions

12-14 nm in diameter Nanodiscs

Storage

Aliquot and store at -20°C or lower.