

Product Information

MemDX™ Recombinant Human CXCR3 Membrane Protein in Virus-Like Particles (MP-VLPs)

Cat. No.: **MPVLP-013**

This product is for research use only and is not intended for diagnostic use.

This product is recombinant Human CXCR3 in VLPs form. This product is produced from mammalian cells by co-expressing the retroviral structural core polyprotein (gag) and the target membrane protein. MP-VLPs display highly-expressed copies of membrane proteins in their native conformation, providing an alternative to membrane protein stable cell lines, membrane preparations, detergent-solubilized proteins and other membrane protein preparation strategies. MP-VLPs can be used for a wide range of applications in antibody production, antibody discovery, antibody characterization, binding assays and functional assays.

Product Specifications

Host Species

Human

Target Protein

CXCR3

Protein Length

Full length

Protein Class

GPCR

TMD

7

Product Description

Application

ELISA; Antibody Production; Antibody Discovery; Antibody Characterization; Binding Assays; Functional Assays

Expression Systems

HEK293 expression system

Protein Format

Membrane Protein-Virus Like Particles (MP-VLPs)

Form

Liquid

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

CXCR3

Full Name

C-X-C motif chemokine receptor 3

Introduction

This gene encodes a G protein-coupled receptor with selectivity for three chemokines, termed CXCL9/Mig (monokine induced by interferon-g), CXCL10/IP10 (interferon-g-inducible 10 kDa protein) and CXCL11/I-TAC (interferon-inducible T cell a-chemoattractant). Binding of chemokines to this protein induces cellular responses that are involved in leukocyte traffic, most notably integrin activation, cytoskeletal changes and chemotactic migration. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the isoforms (CXCR3-B) shows high affinity binding to chemokine, CXCL4/PF4.

Alternative Names

GPR9; MigR; CD182; CD183; Mig-R; CKR-L2; CMKAR3; IP10-R; C-X-C chemokine receptor type 3; G protein-coupled receptor 9; IP-10 receptor; Mig receptor; chemokine (C-X-C motif) receptor 3; chemokine receptor 3; interferon-inducible protein 10 receptor

Gene ID

[2833](#)

UniProt ID

[P49682](#)