

Product Information

MemDX™ Recombinant Human GLP1R Membrane Protein in Virus-Like Particles (MP-VLPs)

Cat. No.: **MPVLP-017**

This product is for research use only and is not intended for diagnostic use.

This product is recombinant Human GLP1R in VLPs form. This product is produced from mammalian cells by co-expressing the retroviral structural core polyprotein (gag) and the target membrane protein. MP-VLPs display highly-expressed copies of membrane proteins in their native conformation, providing an alternative to membrane protein stable cell lines, membrane preparations, detergent-solubilized proteins and other membrane protein preparation strategies. MP-VLPs can be used for a wide range of applications in antibody production, antibody discovery, antibody characterization, binding assays and functional assays.

Product Specifications

Host Species

Human

Target Protein

GLP1R

Protein Length

Full length

Protein Class

GPCR

TMD

7

Product Description

Application

ELISA; Antibody Production; Antibody Discovery; Antibody Characterization; Binding Assays; Functional Assays

Expression Systems

HEK293 expression system

Protein Format

Membrane Protein-Virus Like Particles (MP-VLPs)

Form

Liquid

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

GLP1R

Full Name

Glucagon like peptide 1 receptor

Introduction

This gene encodes a 7-transmembrane protein that functions as a receptor for glucagon-like peptide 1 (GLP-1) hormone, which stimulates glucose-induced insulin secretion. This receptor, which functions at the cell surface, becomes internalized in response to GLP-1 and GLP-1 analogs, and it plays an important role in the signaling cascades leading to insulin secretion. It also displays neuroprotective effects in animal models. Polymorphisms in this gene are associated with diabetes. The protein is an important drug target for the treatment of type 2 diabetes and stroke. Alternative splicing of this gene results in multiple transcript variants.

Alternative Names

GLP-1; GLP-1R; GLP-1-R; glucagon-like peptide 1 receptor; GLP-1 receptor; GLP1 receptor; seven transmembrane helix receptor

Gene ID

[2740](#)

UniProt ID

[P43220](#)