

Product Information

MemDX™ Recombinant Human CD20 Membrane Protein in Virus-Like Particles (MP-VLPs)

Cat. No.: **MPVLP-039**

This product is for research use only and is not intended for diagnostic use.

This product is recombinant Human CD20 in VLPs form. This product is produced from mammalian cells by co-expressing the retroviral structural core polypeptide (gag) and the target membrane protein. MP-VLPs display highly-expressed copies of membrane proteins in their native conformation, providing an alternative to membrane protein stable cell lines, membrane preparations, detergent-solubilized proteins and other membrane protein preparation strategies. MP-VLPs can be used for a wide range of applications in antibody production, antibody discovery, antibody characterization, binding assays and functional assays.

Product Specifications

Host Species

Human

Target Protein

CD20

Protein Length

Full length

Protein Class

B-cell activation

TMD

4

Product Description

Application

ELISA; Antibody Production; Antibody Discovery; Antibody Characterization; Binding Assays; Functional Assays

Expression Systems

HEK293 expression system

Protein Format

Membrane Protein-Virus Like Particles (MP-VLPs)

Form

Liquid

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

CD20

Full Name

Membrane spanning 4-domains A1

Introduction

This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. This gene encodes a B-lymphocyte surface molecule which plays a role in the development and differentiation of B-cells into plasma cells. This family member is localized to 11q12, among a cluster of family members. Alternative splicing of this gene results in two transcript variants which encode the same protein.

Alternative Names

B1; S7; Bp35; CD20; FMC7; CVID5; MS4A2; LEU-16; MS4A1; B-lymphocyte antigen CD20; B-lymphocyte cell-surface antigen B1; CD20 antigen; CD20 receptor; leukocyte surface antigen Leu-16; membrane-spanning 4-domains, subfamily A, member 1

Gene ID

[931](#)

UniProt ID

[P11836](#)