

Product Information

MemDX™ Negative Control Virus-Like Particles (MP-VLPs)

Cat. No.: **MPVLP-050**

This product is for research use only and is not intended for diagnostic use.

This product is negative virus-like particle controls for membrane protein in virus-like particles. This product is produced from HEK293 by expressing the retroviral structural core polypeptide (gag).

Product Specifications

Host Species

Human

Target Protein

VLP

Protein Length

Full length

Product Description

Application

ELISA; Antibody Production; Antibody Discovery; Antibody Characterization; Binding Assays; Functional Assays

Expression Systems

HEK293 expression system

Protein Format

Membrane Protein-Virus Like Particles (MP-VLPs)

Form

Liquid

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

VLP

Alternative Names

VLPs; Virus-like particles