

Product Information

MemDX™ Membrane Protein Pig IL23A (Interleukin 23 subunit alpha) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX0903K**

This product is for research use only and is not intended for diagnostic use.

This product is a Pig IL23A membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Pig

Target Protein

IL23A

Protein Length

Full Length

Protein Class

Immunity

Molecular Weight

26.1 kDa

Sequence

RAVPEGSSPAWAQGQQLSQQQLCTLAHTAHLPMGHVDLPREEGDDETTSEVPHIQCGDGCDPQGLRDNSQSCLQRIHQGLVFYE

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus, Myc tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Reconstitute protein in deionized sterile water

Purity

>85% as determined by SDS-PAGE.

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

IL23A

Full Name

Interleukin 23 subunit alpha

Alternative Names

IL23A; SGRF; il23p19; IL-23 subunit alpha; IL-23-A; IL-23p19; interleukin 23, alpha subunit p19; interleukin-23 p19 subunit; interleukin-23 subunit p19; Interleukin 23 subunit alpha

Gene ID

[100155546](#)

UniProt ID

[Q9N2H9](#)