

Product Information

MemDX™ Membrane Protein Mouse Ccr8 (Chemokine (C-C motif) receptor 8) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX0946K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Ccr8 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Ccr8

Protein Length

Full Length

Protein Class

GPCR

Molecular Weight

42.8 kDa

TMD

7

Sequence

MDYTMEPNVTMTDYYPDFFTAPCDAEFLLRGSMPLYLAILYCVLFVLGLLGNSLVILVLVGCKKLRSITDIYLLNLAASDLLFVLSIPFQT

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Purity

>85% as determined by SDS-PAGE.

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

Ccr8

Full Name

Chemokine (C-C motif) receptor 8

Alternative Names

Ccr8; C-C; CCR-8; CKR-8; Cmkbr; mCCR8; Cmkbr8; CC-CKR-8; C-C CKR-8; C-C motif chemokine receptor 8; Chemokine (C-C motif) receptor 8

Gene ID

[12776](#)

UniProt ID

[P56484](#)