

# Product Information

## MemDX™ Membrane Protein Human UL2 (Protein UL2) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX0992K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human UL2 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

UL2

#### Protein Length

Full Length

#### Protein Class

Virus

#### TMD

1

#### Sequence

MAEDSVAILVEDDNDAYPSFGTLPASHAQYGFRLLRGIFLITLVIWTVVWLKLLRDALL

### Product Description

#### Expression Systems

*in vitro* *E.coli* expression system

#### Tag

10xHis tag at the N-terminus

#### Protein Format

Soluble

#### Form

Liquid or Lyophilized powder

#### Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

### **Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### **Target**

#### **Target Protein**

UL2

#### **Full Name**

Protein UL2

#### **Alternative Names**

UL2; Uncharacterized protein UL2; Protein UL2

#### **Gene ID**

[3077452](#)

#### **UniProt ID**

[P16754](#)