

Product Information

MemDX™ Membrane Protein Human C4orf3 (Chromosome 4 open reading frame 3)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX0993K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human C4orf3 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

C4orf3

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MEVDAPGVDGRDGLRERRGFSEGGQRNQFDVRPQSGANGLPKHSYWLDLWLFILFDVVFVFLFVYFLP

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis and SUMO tag

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

C4orf3

Full Name

Chromosome 4 open reading frame 3

Alternative Names

C4orf3; ALN; HCVFTP1; uncharacterized protein C4orf3; HCV F-transactivated protein 1; another-regulin; hepatitis C virus F protein-transactivated protein 1; Chromosome 4 open reading frame 3

Gene ID

[401152](#)

UniProt ID

[Q8WVX3](#)