

Product Information

MemDX™ Membrane Protein Human SMIM15 (Small integral membrane protein 15)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX0998K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human SMIM15 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SMIM15

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MFDIKAWAEYVVEWAAKDPYGFLT TVILAL TPLFLASAVLSWKLAKMIEAREKEQKKKQKRQENIAKAKRLKKD

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SMIM15

Full Name

Small integral membrane protein 15

Alternative Names

SMIM15; C5orf43; UPF0542 protein C5orf43; Small integral membrane protein 15

Gene ID

[643155](#)

UniProt ID

[Q7Z3B0](#)