

Product Information

MemDX™ Membrane Protein Mouse Tmem234 (Transmembrane protein 234) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1146K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Tmem234 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Tmem234

Protein Length

Full Length

Protein Class

Receptor

TMD

3

Sequence

MAASWGQVLALVLVAALWGGTQPLLKRASSGLEQVRERTWAWQLLQEIKALFGNTEYLMPFLLNQSGSLLYYLTLASTDLTLAVPI

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

Tmem234

Full Name

Transmembrane protein 234

Alternative Names

Tmem234; 1500002D11Rik; 2510006D16Rik; 4933407D05Rik; Transmembrane protein 234

Gene ID

[76799](#)

UniProt ID

[Q8R1E7](#)