

# Product Information

## **MemDX™ Membrane Protein Mouse Tmem242 (Transmembrane protein 242) Expressed *in vitro* E.coli expression system, Full Length**

Cat. No.: **MPX1147K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Tmem242 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Mouse

#### Target Protein

Tmem242

#### Protein Length

Full Length

#### Protein Class

Receptor

#### TMD

2

#### Sequence

METSGPGPGESSELEAPGSPDDRFLVKGGIFLGSAAGMLAGFVTTLAKKKSPWFNKGTMATAALPESGSSLALRALGWG

### Product Description

#### Expression Systems

*in vitro* E.coli expression system

#### Tag

10xHis tag at the N-terminus

#### Protein Format

Soluble

#### Form

Liquid or Lyophilized powder

#### Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

### **Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### **Target**

#### **Target Protein**

Tmem242

#### **Full Name**

Transmembrane protein 242

#### **Alternative Names**

Tmem242; 1110008A10Rik; 2310046K16Rik; 5730437N04Rik; Transmembrane protein 242

#### **Gene ID**

[70544](#)

#### **UniProt ID**

[Q8VCR3](#)