

Product Information

MemDX™ Membrane Protein Human OPALIN (Oligodendrocytic myelin paranodal and inner loop protein) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1152K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human OPALIN membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

OPALIN

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MSFSLNFTLPANTTSSPVTGGKETDCGPSLGLAAGIPLL VATALLVALLFTLIHRRRSSIEAMEESDRPCEISEIDDNPKISENPRRSPT

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

OPALIN

Full Name

Oligodendrocytic myelin paranodal and inner loop protein

Alternative Names

OPALIN; TMP10.htmP10; TMEM10; opalin; transmembrane protein 10; transmembrane protein TMP10; Oligodendrocytic myelin paranodal and inner loop protein

Gene ID

[93377](#)

UniProt ID

[Q96PE5](#)