

# Product Information

## **MemDX™ Membrane Protein Mouse Ifitm2 (Interferon induced transmembrane protein 2)**

### **Expressed *in vitro* *E.coli* expression system, Full Length**

Cat. No.: **MPX1167K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Ifitm2 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### **Product Specifications**

#### **Host Species**

Mouse

#### **Target Protein**

Ifitm2

#### **Protein Length**

Full Length

#### **Protein Class**

Immunity

#### **TMD**

2

#### **Sequence**

MSHNSQAFLSTNAGLPPSYETIKEEYGVTELGEPSNSAVVRTTVINMPREVSVPDHVVWSLFNTLFFNACCLGFVAYAYSVKSRDR

### **Product Description**

#### **Expression Systems**

*in vitro* *E.coli* expression system

#### **Tag**

10xHis tag at the N-terminus

#### **Protein Format**

Soluble

#### **Form**

Liquid or Lyophilized powder

#### **Buffer**

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

### **Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### **Target**

#### **Target Protein**

Ifitm2

#### **Full Name**

Interferon induced transmembrane protein 2

#### **Alternative Names**

Ifitm2; fra; mil; Ifit; mil-3; DSPA2c; Ifitm3l; fragilis3; dispanin subfamily A member 2c; fragilis R; fragilis protein 3; interferon induced transmembrane protein, like; Interferon induced transmembrane protein 2

#### **Gene ID**

[80876](#)

#### **UniProt ID**

[Q99J93](#)