

Product Information

MemDX™ Membrane Protein Human PLGRKT (Plasminogen receptor with a C-terminal lysine) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1172K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human PLGRKT membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PLGRKT

Protein Length

Full Length

Protein Class

Receptor

TMD

2

Sequence

MGFIFSKSMNESMKNQKEFMLMNARLQLERQLIMQSEMRERQMAMQIAWSREFLYFGTFFGLAAISLTAGAIKKKKPAFLVPIVPL

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

PLGRKT

Full Name

Plasminogen receptor with a C-terminal lysine

Alternative Names

PLGRKT; AD025; MDS030; C9orf46; PLG-RKT; Plg-R(KT); plasminogen receptor (KT); 5033414D02Rik; plasminogen receptor, C-terminal lysine transmembrane protein; transmembrane protein C9orf46; Plasminogen receptor with a C-terminal lysine

Gene ID

[55848](#)

UniProt ID

[Q9HBL7](#)