

Product Information

MemDX™ Membrane Protein Human SHISAL2A (Shisa like 2A) Expressed *in vitro* E.coli **expression system, Full Length**

Cat. No.: **MPX1310K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human SHISAL2A membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SHISAL2A

Protein Length

Full Length

Protein Class

Receptor

TMD

2

Sequence

MSGACTSYVSAEQEVVRGFSCPRPGGEAAVFCCGFRDHKYCCDDPHSFFPYEHSYMWLWSIGALIGLSVAADVLLAFIVTACVLC

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

SHISAL2A

Full Name

Shisa like 2A

Alternative Names

SHISAL2A; FAM159A; PRO7171; WWLS2783; protein shisa-like-2A; family with sequence similarity 159 member A; membrane protein FAM159A; Shisa like 2A

Gene ID

[348378](#)

UniProt ID

[Q6UWV7](#)