

Product Information

MemDX™ Membrane Protein Mouse Fam210b (Family with sequence similarity 210, member B) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1315K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Fam210b membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Fam210b

Protein Length

Full Length

Protein Class

Receptor

TMD

2

Sequence

MAGLLTLLGPAGRVSTRRLRPLAPWLLGTATSCAPPLWALALSHVPDARLLRTARGDCLSRQEPNRTPEPGGGVTGTEKKLSRTQ

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Fam210b

Full Name

Family with sequence similarity 210, member B

Alternative Names

Fam210b; AI430961; AV108736; 2010011I20Rik; protein FAM210B, mitochondrial; transmembrane protein C20orf108 homolog; Family with sequence similarity 210, member B

Gene ID

[67017](#)

UniProt ID

[Q9D8B6](#)