

Product Information

MemDX™ Membrane Protein Mouse Rnf185 (Ring finger protein 185) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1318K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Rnf185 membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Rnf185

Protein Length

Full Length

Protein Class

Transferase

TMD

2

Sequence

MASKGPSASASTENSNAGGPSGSSNGTGESGGQDSTFECNICLDTAKDAVISLCGHLFCWPCLHQWLETRPNRQVCPVCKAGISF

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

Rnf185

Full Name

Ring finger protein 185

Alternative Names

Rnf185; E3 ubiquitin-protein ligase RNF185; RING-type E3 ubiquitin transferase RNF185; Ring finger protein 185

Gene ID

[193670](#)

UniProt ID

[Q91YT2](#)