

Product Information

MemDX™ Membrane Protein Mouse Mpv17I2 (MPV17 mitochondrial membrane protein-like

2) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1340K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Mpv17I2 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Mpv17I2

Protein Length

Full Length

Protein Class

Receptor

TMD

3

Sequence

MALGGWRWARKALAAGRPLFQGRALLLTNTLGCGVLMAGDGARQVWEVRARPGQRFSARRSASMFAVGCSMGPFLLHFWYLV

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Mpv17l2

Full Name

MPV17 mitochondrial membrane protein-like 2

Alternative Names

Mpv17l2; Fksg2; Fksg24; mpv17-like protein 2; MPV17 mitochondrial membrane protein-like 2

Gene ID

[234384](#)

UniProt ID

[Q8VIK2](#)