

Product Information

MemDX™ Membrane Protein Mouse Rnf152 (Ring finger protein 152) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1347K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Rnf152 membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Rnf152

Protein Length

Full Length

Protein Class

Transferase

TMD

1

Sequence

METLSQDSLLEQCICFNYYSPRRRPKLLDCKHTCCSVCLQQMRTSQKDVRCPWCRGITKLPPGFSVSQLPDDPEVLAVIAIPHTSEH

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Rnf152

Full Name

Ring finger protein 152

Alternative Names

Rnf152; A930029B02Rik; E3 ubiquitin-protein ligase RNF152; RING-type E3 ubiquitin transferase RNF152; Ring finger protein 152

Gene ID

[320311](#)

UniProt ID

[Q8BG47](#)