

Product Information

MemDX™ Membrane Protein Human MPZL3 (Myelin protein zero like 3) Expressed *in vitro* *E.coli* expression system, Full Length of Mature Protein

Cat. No.: **MPX1350K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human MPZL3 membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MPZL3

Protein Length

Full Length of Mature Protein

Protein Class

Cell adhesion

TMD

1

Sequence

LEIRADAHVVRGYVGEKIKLKCTFKSTSDVTDKLTIDWTYRPPSSSHTVSIFHYQSFQYPTTAGTFRDRISWVGNNVYKGDASISISNPTI

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

MPZL3

Full Name

Myelin protein zero like 3

Alternative Names

MPZL3; myelin protein zero-like protein 3; Myelin protein zero like 3

Gene ID

[196264](#)

UniProt ID

[Q6UWV2](#)