

Product Information

MemDX™ Membrane Protein Human VSTM1 (V-set and transmembrane domain containing

1) Expressed *in vitro* E.coli expression system, Full Length of Mature Protein

Cat. No.: **MPX1397K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human VSTM1 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

VSTM1

Protein Length

Full Length of Mature Protein

Protein Class

Cytokine

TMD

1

Sequence

YEDEKKNEKPPKPSLHAWPSSVVEAESNVTLKCQAHSQNVTFVLRKVNDSGYKQEQSSAENAEFPFTDLKPKDAGRYFCAYKTT

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

VSTM1

Full Name

V-set and transmembrane domain containing 1

Alternative Names

VSTM1; SIRL1; SIRL-1; UNQ3033; V-set and transmembrane domain-containing protein 1; LAIR homolog; OSCAR-like transcript-1; signal inhibitory receptor on leukocytes-1; V-set and transmembrane domain containing 1

Gene ID

[284415](#)

UniProt ID

[Q6UX27](#)